

Fluoromount

Cat. No. GTX03629

Application IHC-P

Package
25 ml

PRODUCT

Summary

Fluoromount is an aqueous mounting medium for cover slipping tissues stained with various fluorescent dyes. This product is compatible with FITC, phycoerythrin, phycocyanin and allophycocyanin. It can be used with other fluorescent markers such as Texas Red and rhodamine.

PROPERTIES

Form Liquid

Storage Store at room temperature.

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).