

# Fluoromount/Plus

**Cat. No. GTX03630**

**Application** IHC-P

**Package**  
25 ml

## PRODUCT

### Summary

Fluoromount/Plus is an aqueous-based mounting medium for applying coverslips to tissues stained with fluorescent dyes. This product is compatible with FITC, phycoerythrin, phycocyanin and allophycocyanin. It can also be used with other fluorescent markers such as Texas Red and rhodamine.

## PROPERTIES

**Form** Liquid

**Storage** Store at room temperature.

**Conjugation** Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).