

# Tween 20

**Cat. No. GTX03634**

**Package**  
125 ml

## PRODUCT

### Summary

Tween 20 is a detergent (surfactant) commonly added to buffers and reagents for immunohistochemistry. Its purpose is to decrease background staining and enhance reagent spreading in automated and manual procedures. It is suggested to use the product at concentrations between 0.05% and 0.1% in buffers and reagents.

## PROPERTIES

**Form** Liquid

**Storage** Store at room temperature.

**Conjugation** Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).