

Mouse Tissue Array III - Normal Tissues

Cat. No. GTX03659

Applications IHC-Fr, ISH**Species** Mouse**Package**

5 slide

PRODUCT

Summary

Tissues were snap frozen in liquid nitrogen immediately after excision, and embedded in OCT. All Genetex's frozen tissue products, including frozen tissue arrays, frozen tissue sections, and frozen tissue panels, can be used for both immunohistochemistry and in situ hybridization. These products are ideal for rapid cellular localization of RNA and protein expression, especially for antibodies that don't recognize paraffin embedded tissue. At least one of the tissue slides from each lot was stained with H & E to ensure the quality. Frozen Tissue Section Panel contains 5 to 8 different tissues on one slide. The tissue panels are 5-10 μ m in thickness, mounted on positively charged glass slides, and fixed by cold acetone. Note: Each individual tissue spot should be >90% core retention and total core lost from 1 slide should be no more than 5%.

Tissue type Normal

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

IHC-Fr

Assay dependent

ISH

Assay dependent

Not tested in other applications.

Properties

Storage

Half a year from the date of receipt under proper storage condition

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).