

Mouse anti-Human IgG4 pFc' antibody [HP6023] (HRP)

Cat. No. GTX03663

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Applications	WB, IHC-Fr, FCM, IP, ELISA
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	1:2000-1:8000

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 50% glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG4 myeloma protein
Purification	Purified IgG3
Conjugation	Horseradish peroxidase(HRP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX03663 ELISA Image

ELISA analysis of purified human IgG1, IgG2, IgG3, and IgG4 using serially diluted GTX03663 Mouse anti-Human IgG4 pFc' antibody [HP6023] (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).