

Gasdermin C antibody

Cat. No. GTX03667

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG

Package
100 µl

Applications

Calculated MW 58 kDa. ([Note](#))

Product Note We do not recommend use of this product for Sheep,Bovine samples.

Properties

Form Liquid

Buffer PBS, 150mM NaCl, 50% Glycerol

Preservative 0.02% sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Immunogen A synthesized peptide derived from human GSDMC(Accession Q9BYG8), corresponding to amino acid residues Q120-F144.

Purification Purified by antigen-affinity chromatography
From serum

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).