

Goat anti-Rabbit IgA antibody (HRP)

Cat. No. GTX03700-01

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Rabbit

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	1:10000-1:100000

Not tested in other applications.

Product Note

This antibody recognizes rabbit IgA and shows no cross-reactivity with other rabbit immunoglobulin classes in immunoelectrophoresis.

Properties

Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	0.05% proclin400
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified rabbit IgA.
Purification	Purified by affinity chromatography From serum
Conjugation	Horseradish peroxidase(HRP)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).