

Mouse Heart (Normal) tissue lysate

Cat. No. GTX03710**Species** Mouse**Package**

100 µg

PRODUCT

Summary

Tissue total protein is prepared from whole tissue homogenates. It presents a consistent pattern on SDS-PAGE analysis. The total protein is stored in a buffer with protease inhibitor cocktail.

Tissue type

Normal

Properties

Form

Liquid

BufferHEPES, MgCl₂, KCl, EDTA, Sucrose, Glycerol, Sodium deoxycholate, NP-40, and a cocktail of protease inhibitors**Preservative**

No preservatives

Storage

Store the protein lysate at -70°C for long term storage. After the protein is thawed, store it at 2 - 8°C for continuous use. Avoid freeze and thaw cycles.

Concentration

5 mg/ml (Please refer to the vial label for the specific concentration.)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).