

Human PKC iota protein, His tag

Cat. No. GTX03726-pro
Species Human

Package
50 µg

Properties
Form Liquid

Buffer PBS, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence 1-596aa

Expression System HEK293 cells


Purification >80% by SDS-PAGE

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX03726-pro Image

SDS-PAGE of GTX03726-pro Human PKC iota protein, His tag.



For full product information, images and publications, please visit our [website](#).