

Lymphocytic Choriomeningitis virus (LCMV) Nucleoprotein antibody [VL-4]

Cat. No. GTX03737

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	ICC/IF, FCM
Reactivity	Lymphocytic Choriomeningitis Virus

Package

500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Product Note

This antibody has been shown to stain LCMV-infected cell internally with no surface staining. This antibody does not react with vaccinia, vesicular stomatitis or influenza virus-infected cells in the case of internal or surface staining.

Properties

Form	Liquid
Buffer	Filter-sterilized PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	LCMV strain WE
Purification	Protein G purified From tissue culture supernatant
Purity	>95% by SDS-PAGE
Endotoxin	< 0.002 EU/µg (determined by LAL assay)
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 16 Page 1 of 2

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).