

Human IgM isotype control

Cat. No. GTX03790

Host	Human
Clonality	Polyclonal
Isotype	IgM
Applications	ELISA

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	$\leq 2\mu\text{g}$

Not tested in other applications.

Properties

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified IgM From serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx03790 ELISA Image

ELISA analysis of serially diluted GTx03790 Human IgM isotype control using Mouse Anti-Human IgG Hinge, Mouse Anti-Human IgG2 Fc, Mouse Anti-Human IgG3 Hinge, Mouse Anti-Human IgG4 PFc', Mouse Anti-Human IgA1, Mouse Anti-Human IgA2, Mouse Anti-Human IgD, Mouse Anti-Human IgE Fc, Mouse Anti-Human IgM, Mouse Anti-Human Kappa, and Mouse Anti-Human Lambda followed by Streptavidin-HRP and quantified.



For full product information, images and publications, please visit our [website](#).