

Human IgG1 Fc

Cat. No. GTX03801-pro

Species Human

Package
500 µg

Properties

Form Liquid

Buffer Filter-sterilized PBS

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence Recombinant Human IgG1 Fc

Expression System Purified from NS0 cell medium

Purification Protein A purified
From tissue culture supernatant

Purity >95% by SDS-PAGE

Endotoxin <2.0EU/1mg (determined by the LAL method)

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).