

Human IgG1 Fc

Cat. No. GTX03801-pro

Species	Human
----------------	-------

Package 500 µg

Properties

Form	Liquid
-------------	--------

Buffer	Filter-sterilized PBS
---------------	-----------------------

Preservative	No preservatives
---------------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
----------------	--

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
----------------------	--

Region/Sequence	Recombinant Human IgG1 Fc
------------------------	---------------------------

Expression System	Purified from NS0 cell medium
--------------------------	-------------------------------

Purification	Protein A purified From tissue culture supernatant
---------------------	---

Purity	>95% by SDS-PAGE
---------------	------------------

Endotoxin	<2.0EU/1mg (determined by the LAL method)
------------------	---

Conjugation	Unconjugated
--------------------	--------------

Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>
-------------	---



For full product information, images and publications, please visit our [website](#).