

Human NT-proBNP1 protein

Cat. No. GTX04089-pro

Species Human

Package

50 µg

Properties

Form Lyophilized powder

Buffer Reconstitute the lyophilized protein in sterile distilled H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions. Lyophilized from a 0.2µm filtered concentrated solution in 20mM Tris-HCl, pH 8.0 and 150mM NaCl.

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Region/Sequence HPLGSPGSAS DLETSGLQEQ RNHLQGKLSE LQVEQTSLEP LQESPRPTGV WKSREVATEG IRGHRKMVLY TLRAPR

Expression System E. coli

Purification Purified by proprietary chromatographic techniques

Purity > 98% by SDS-PAGE and RP-HPLC

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 1