

Mouse Anti-Armenian Hamster IgG2+IgG3 antibody [SB139e] (Biotin)

Cat. No. GTX04118-02

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM, ELISA
Reactivity	Armenian Hamster

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
ELISA	1:5000-1:10000

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Armenian Hamster IgG
Purification	Purified IgG1
Conjugation	Biotin


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04118-02 ELISA Image

ELISA analysis of purified Armenian Hamster immunoglobulins using serially diluted GTX04118-02 Mouse Anti-Armenian Hamster IgG2+IgG3 antibody [SB139e] (Biotin) followed by Streptavidin (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).