

Mouse Anti-Armenian Hamster IgG2+IgG3 antibody [SB139e] (PE)

Cat. No. GTX04118-08

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM, ELISA
Reactivity	Armenian Hamster

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	$\leq 0.1 \mu\text{g}/10^6\text{cells}$
ELISA	$\leq 1 \mu\text{g}/\text{mL}$

Note : The suggested use of these reagents is in a final volume of 100µl.

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS with a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Armenian Hamster IgG
Purification	Purified IgG1
Conjugation	Phycoerythrin (PE) Wavelength

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx04118-08 ELISA Image

ELISA analysis of purified Armenian Hamster immunoglobulins using serially diluted GTx04118-08 Mouse Anti-Armenian Hamster IgG2+IgG3 antibody [SB139e] (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).