

Goat Anti-Human IgG (Fc) antibody, pre-adsorbed (PE)

Cat. No. GTX04125-08

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	≤ 1 µg/mL

Not tested in other applications.

Product Note Pre-adsorbed with Human IgG Fab, IgM, and IgA. May react with IgG from other species.

Properties

Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from goats hyperimmunized with human IgG
Purification	Purified by antigen-affinity chromatography
Conjugation	Phycoerythrin (PE) Wavelength


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04125-08 ELISA Image

ELISA analysis of human immunoglobulins using serially diluted GTx04125-08 Goat Anti-Human IgG (Fc) antibody, pre-adsorbed (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).