

Goat Anti-Human lambda light chain antibody, pre-adsorbed (AP)

Cat. No. GTX04134-03

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	FCM, ELISA
Reactivity	Human

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
ELISA	1:2000-1:4000

Not tested in other applications.

Product Note

Pre-adsorbed with Human κ light chains; mouse immunoglobulins and pooled sera. May react with λ light chains from other species.

Properties

Form	Liquid
Buffer	50mM Tris, 1mM MgCl ₂ , 50% Glycerol
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from goats hyperimmunized with human λ light chains
Purification	Purified by antigen-affinity chromatography
Conjugation	Alkaline Phosphatase (AP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04134-03 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTX04134-03 Goat Anti-Human lambda light chain antibody, pre-adsorbed (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).