

Goat Anti-Rat IgM (Mu chain) antibody, pre-adsorbed (PE)

Cat. No. GTX04137-08

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Application	WB, ICC/IF, IHC-P, IHC-Fr, FACS, ELISA
Reactivity	Rat

Package
500 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FACS	Assay dependent
ELISA	≤ 1 µg/mL

Not tested in other applications.

Product Note Pre-adsorbed with Rat IgG. May react with immunoglobulins from other species.

PROPERTIES

Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from goats hyperimmunized with rat IgM
Purification	Purified by antigen-affinity chromatography
Conjugation	Phycoerythrin (PE)




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx04137-08 ELISA Image

ELISA analysis of purified rat immunoglobulins using serially diluted GTx04137-08 Goat Anti-Rat IgM (Mu chain) antibody, pre-adsorbed (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).