

Mouse Anti-Rat IgG2b antibody [2B10A8] (AP)

Cat. No. GTX04141-03

| | |
|--------------|------------------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG2b |
| Applications | WB, IHC-Fr, FCM, ELISA |
| Reactivity | Rat |

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB | Assay dependent |
| IHC-Fr | Assay dependent |
| FCM | Assay dependent |
| ELISA | 1:1000-1:4000 |

Not tested in other applications.

Properties

| | |
|---------------|---|
| Form | Liquid |
| Buffer | 50mM Tris, 1mM MgCl ₂ , 50% Glycerol |
| Preservative | 0.09% sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.) |
| Immunogen | Rat IgG2b hybridoma |
| Purification | Purified IgG2b |
| Conjugation | Alkaline Phosphatase (AP) |


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04141-03 ELISA Image

ELISA analysis of purified rat immunoglobulins using serially diluted GTX04141-03 Mouse Anti-Rat IgG2b antibody [2B10A8] (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).