

# Goat Anti-Rabbit IgM (Mu chain) antibody, pre-adsorbed (PE)

**Cat. No. GTX04146-08**

<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, IHC-P, IHC-Fr, FCM, ELISA, ELISPOT
<b>Reactivity</b>	Rabbit

**Package**  
500 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
ELISA	≤ 1 µg/mL
ELISPOT	Assay dependent

Not tested in other applications.

**Product Note** Pre-adsorbed with Rabbit IgG. May react with immunoglobulins from other species.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, a stabilizer
<b>Preservative</b>	0.09% sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
<b>Concentration</b>	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Pooled antisera from goats hyperimmunized with rabbit IgM
<b>Purification</b>	Purified by antigen-affinity chromatography
<b>Conjugation</b>	Phycoerythrin (PE) <a href="#">Wavelength</a>




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

**DATA IMAGES**

**GTX04146-08 ELISA Image**

ELISA analysis of purified rabbit immunoglobulins using serially diluted GTX04146-08 Goat Anti-Rabbit IgM (Mu chain) antibody, pre-adsorbed (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).