

Goat Anti-Rabbit IgG (Fab) antibody, pre-adsorbed (HRP)

Cat. No. GTX04147-01

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Rabbit

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	1:2000-1:4000

Not tested in other applications.

Product Note Pre-adsorbed with Rabbit IgG Fc. May react with immunoglobulins from other species.

Properties


Form	Liquid
Buffer	PBS, 50% glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from goats hyperimmunized with rabbit IgG
Purification	Purified by antigen-affinity chromatography
Conjugation	Horseradish peroxidase(HRP)

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04147-01 ELISA Image

ELISA analysis of purified rabbit immunoglobulins using serially diluted GTx04147-01 Goat Anti-Rabbit IgG (Fab) antibody, pre-adsorbed (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).