

Mouse Anti-Rabbit (Light chain) antibody [SB62a]

Cat. No. GTX04152

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ELISA
Reactivity	Rabbit

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	$\leq 1 \mu\text{g/mL}$

Not tested in other applications.

Properties

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit light chains
Purification	Purified IgG2b
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx04152 ELISA Image

ELISA analysis of purified rabbit immunoglobulins using serially diluted GTx04152 Mouse Anti-Rabbit (Light chain) antibody [SB62a] followed by Goat Anti-Mouse IgG2b antibody, pre-absorbed (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).