

Mouse Anti-Rabbit (Light chain) antibody [SB62a] (AP)

Cat. No. GTX04152-03

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ELISA
Reactivity	Rabbit

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	1:1000-1:2000

Not tested in other applications.

Properties

Form	Liquid
Buffer	50mM Tris, 1mM MgCl ₂ , 50% Glycerol
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit light chains
Purification	Purified IgG2b
Conjugation	Alkaline Phosphatase (AP)

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx04152-03 ELISA Image

ELISA analysis of purified rabbit immunoglobulins using serially diluted GTx04152-03 Mouse Anti-Rabbit (Light chain) antibody [SB62a] (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).