

Mouse Anti-Rabbit IgG (Fc) antibody [2A9] (PE)

Cat. No. GTX04153-08

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, ELISA
Reactivity	Rabbit

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
FCM	$\leq 0.3 \mu\text{g}/10^6\text{cells}$
ELISA	$\leq 1 \mu\text{g/mL}$

Note : The suggested use of these reagents is in a final volume of 100µl.

Not tested in other applications.

Product Note Minimal reactivity to human, mouse, rat, goat, bovine, horse, donkey, sheep, guinea pig, and chicken serum.

Properties


Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant rabbit IgG Fc fragment
Purification	Purified IgG1
Conjugation	Phycoerythrin (PE) Wavelength

For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX04153-08 ELISA Image

ELISA analysis of purified rabbit immunoglobulins using serially diluted GTX04153-08 Mouse Anti-Rabbit IgG (Fc) antibody [2A9] (PE).


GTX04153-08 FCM Image

FACS analysis of BALB/c mouse splenocytes using Rabbit Anti-Mouse IgM antibody followed by GTX04153-08 Mouse Anti-Rabbit IgG (Fc) antibody [2A9] (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).