

Rabbit Anti-Human IgG antibody, F(ab')₂ fragment, pre-adsorbed (AP)

Cat. No. GTX04154-03

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG F(ab') ₂
Applications	WB, ICC/IF, ELISA
Reactivity	Human

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
ELISA	1:2000-1:4000

Not tested in other applications.

Product Note

Reacts with the heavy and light chains of human IgG and the light chains of human IgM and IgA. Pre-adsorbed with None. May react with immunoglobulins from other species

Properties

Form	Liquid
Buffer	50mM Tris, 1mM MgCl ₂ , 50% Glycerol
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Pepsin digest of Rabbit Anti-Human IgG(H+L).
Purification	Purified IgG
Conjugation	Alkaline Phosphatase (AP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04154-03 ELISA Image

ELISA analysis of purified human immunoglobulins using GTX04154-03 Rabbit Anti-Human IgG antibody, F(ab')₂ fragment, pre-adsorbed (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).