

# Rabbit Anti-Sheep IgG antibody, F(ab')<sub>2</sub> fragment, pre-adsorbed

**Cat. No. GTX04157**

|                     |                         |
|---------------------|-------------------------|
| <b>Host</b>         | Rabbit                  |
| <b>Clonality</b>    | Polyclonal              |
| <b>Isotype</b>      | IgG F(ab') <sub>2</sub> |
| <b>Applications</b> | WB, ELISA               |
| <b>Reactivity</b>   | Sheep                   |

**Package**  
500 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | Assay dependent      |
| ELISA              | Assay dependent      |

Not tested in other applications.

### Product Note

Pre-adsorbed with Human serum proteins (SP). May react with immunoglobulins from other species and the light chains of other sheep immunoglobulins.

## Properties

|                      |                                                                                         |
|----------------------|-----------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                  |
| <b>Buffer</b>        | Borate buffered saline                                                                  |
| <b>Preservative</b>  | No preservatives                                                                        |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| <b>Concentration</b> | 0.5 mg/ml (Please refer to the vial label for the specific concentration.)              |
| <b>Immunogen</b>     | Pepsin digest of Rabbit Anti-Sheep IgG(H+L).                                            |
| <b>Purification</b>  | Purified IgG                                                                            |
| <b>Conjugation</b>   | Unconjugated                                                                            |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx04157 ELISA Image

ELISA analysis of purified sheep IgG and human serum using GTx04157 Rabbit Anti-Sheep IgG antibody, F(ab')<sub>2</sub> fragment, pre-adsorbed followed by Goat Anti-Rabbit IgG antibody, pre-absorbed (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).