

Rabbit Anti-Goat IgG (Fab) antibody, pre-adsorbed (Biotin)

Cat. No. GTX04159-02

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Goat

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	1:5000-1:20000

Not tested in other applications.

Product Note

Pre-adsorbed with Goat IgG Fc. May react with immunoglobulins from other species and the light chains of other goat immunoglobulins.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from rabbits hyperimmunized with goat IgG
Purification	Purified by antigen-affinity chromatography
Conjugation	Biotin

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 13 Page 1 of 2

DATA IMAGES

GTX04159-02 ELISA Image

ELISA analysis of purified goat immunoglobulins using GTX04159-02 Rabbit Anti-Goat IgG (Fab) antibody, pre-adsorbed (Biotin) followed by Streptavidin (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).