

## Goat Anti-Golden Syrian Hamster IgG antibody, pre-adsorbed (PE)

Cat. No. GTX04160-08

|              |                               |
|--------------|-------------------------------|
| Host         | Goat                          |
| Clonality    | Polyclonal                    |
| Isotype      | IgG                           |
| Applications | WB, IHC-P, IHC-Fr, FCM, ELISA |
| Reactivity   | Golden Syrian Hamster         |

Package  
500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution                    |
|--------------------|-----------------------------------------|
| WB                 | Assay dependent                         |
| IHC-P              | Assay dependent                         |
| IHC-Fr             | Assay dependent                         |
| FCM                | $\leq 0.1 \mu\text{g}/10^6\text{cells}$ |
| ELISA              | $\leq 1 \mu\text{g/mL}$                 |

**Note :** The suggested use of these reagents is in a final volume of 100µl.

Not tested in other applications.

## Product Note

Pre-adsorbed with Mouse and rat immunoglobulins and pooled sera. May react with immunoglobulins from other species and the light chains of other hamster immunoglobulins.

## Properties


|               |                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                     |
| Buffer        | PBS, a stabilizer                                                                                                          |
| Preservative  | 0.09% sodium azide                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light. |
| Concentration | 0.5 mg/ml (Please refer to the vial label for the specific concentration.)                                                 |
| Immunogen     | Pooled antisera from goats hyperimmunized with hamster IgG                                                                 |
| Purification  | Purified by antigen-affinity chromatography                                                                                |
| Conjugation   | Phycoerythrin (PE) <a href="#">Wavelength</a>                                                                              |

For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

**DATA IMAGES**

**GTX04160-08 FCM Image**

FACS analysis of BALB/c mouse splenocytes using Golden Syrian Hamster Anti-Mouse CD3 epsilon antibody followed by GTX04160-08 Goat Anti-Golden Syrian Hamster IgG antibody, pre-adsorbed (PE).


**GTX04160-08 ELISA Image**

ELISA analysis of purified Golden Syrian Hamster IgG using GTX04160-08 Goat Anti-Golden Syrian Hamster IgG antibody, pre-adsorbed (PE).


**GTX04160-08 FCM Image**

FACS analysis of BALB/c mouse thymocytes using Golden Syrian Hamster Anti-Mouse CD28 antibody followed by GTX04160-08 Goat Anti-Golden Syrian Hamster IgG antibody, pre-adsorbed (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).