

Goat Anti-Golden Syrian Hamster IgG antibody, F(ab')₂ fragment, pre-adsorbed**Cat. No. GTX04161**

Host	Goat
Clonality	Polyclonal
Isotype	IgG F(ab') ₂
Applications	WB, FCM, ELISA, IHC
Reactivity	Golden Syrian Hamster

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	Assay dependent
ELISA	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Product Note

Pre-adsorbed with Mouse and rat immunoglobulins and pooled sera. May react with immunoglobulins from other species and the light chains of other hamster immunoglobulins.

Properties

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pepsin digest of Goat Anti-Hamster IgG(H+L).
Purification	Purified IgG
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04161 ELISA Image

ELISA analysis of purified Golden Syrian Hamster, mouse, rat immunoglobulins using GTX04161 Goat Anti-Golden Syrian Hamster IgG antibody, F(ab')₂ fragment, pre-adsorbed followed by Swine Anti-Goat IgG(H+L) antibody, pre-absorbed (HRP).


GTX04161 FCM Image

FACS analysis of BALB/c mouse splenocytes using Golden Syrian Hamster Anti-Mouse CD3 epsilon antibody followed by GTX04161 Goat Anti-Golden Syrian Hamster IgG antibody, F(ab')₂ fragment, pre-adsorbed and Swine Anti-Goat IgG antibody, pre-absorbed (FITC).



For full product information, images and publications, please visit our [website](https://www.genetex.com).