

# Goat Anti-Golden Syrian Hamster IgG antibody, F(ab')<sub>2</sub> fragment, pre-adsorbed (PE)

**Cat. No. GTX04161-08**

<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG F(ab') <sub>2</sub>
<b>Applications</b>	FCM, ELISA, IHC
<b>Reactivity</b>	Golden Syrian Hamster

**Package**

250 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	≤ 0.2 µg/10 <sup>6</sup> cells
ELISA	≤ 1 µg/mL
IHC	Assay dependent

**Note : The suggested use of these reagents is in a final volume of 100µl.**

Not tested in other applications.

### Product Note

Pre-adsorbed with Mouse and rat immunoglobulins and pooled sera. May react with immunoglobulins from other species and the light chains of other hamster immunoglobulins.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, a stabilizer
<b>Preservative</b>	0.09% sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
<b>Concentration</b>	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Pepsin digest of Goat Anti-Hamster IgG(H+L).
<b>Purification</b>	Purified IgG
<b>Conjugation</b>	Phycoerythrin (PE) <a href="#">Wavelength</a>




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

**DATA IMAGES**

**GTx04161-08 ELISA Image**

ELISA analysis of purified Golden Syrian Hamster IgG using GTx04161-08 Goat Anti-Golden Syrian Hamster IgG antibody, F(ab')<sub>2</sub> fragment, pre-adsorbed (PE).


**GTx04161-08 FCM Image**

FACS analysis of BALB/c mouse splenocytes using Golden Syrian Hamster Anti-Mouse CD3 epsilon antibody followed by GTx04161-08 Goat Anti-Golden Syrian Hamster IgG antibody, F(ab')<sub>2</sub> fragment, pre-adsorbed (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).