

Goat Anti-Chicken IgY antibody, pre-adsorbed (PE)

Cat. No. GTX04162-08

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Chicken

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
FCM	Assay dependent
ELISA	≤ 1 µg/mL

Not tested in other applications.

Product Note

Pre-adsorbed with None. May react with immunoglobulins from other species and the light chains of other chicken immunoglobulins.

Properties

Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from goats hyperimmunized with IgY
Purification	Purified by antigen-affinity chromatography
Conjugation	Phycoerythrin (PE) Wavelength




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX04162-08 ELISA Image

ELISA analysis of purified chicken immunoglobulins using GTX04162-08 Goat Anti-Chicken IgY antibody, pre-adsorbed (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).