

Rabbit Anti-Mouse IgG antibody, F(ab')₂ fragment, pre-adsorbed

Cat. No. GTX04163

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG F(ab') ₂
Applications	WB, ICC/IF, FCM, ELISA
Reactivity	Mouse

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Product Note

Reacts with the heavy and light chains of mouse IgG and the light chains of mouse IgM and IgA. Pre-adsorbed with None. May react with immunoglobulins from other species

Properties

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pepsin digest of Rabbit Anti-Mouse IgG(H+L).
Purification	Purified IgG
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04163 ELISA Image

ELISA analysis of purified mouse immunoglobulins using GTX04163 Rabbit Anti-Mouse IgG antibody, F(ab')₂ fragment, pre-adsorbed followed by Goat Anti-Rabbit IgG antibody, pre-adsorbed (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).