

Rabbit Anti-Rat IgG antibody, F(ab')₂ fragment, pre-adsorbed (TRITC)

Cat. No. GTX04164-25

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG F(ab') ₂
Applications	WB, ICC/IF, IHC-P, IHC-Fr, ELISA
Reactivity	Rat

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
ELISA	1:100-1:200

Not tested in other applications.

Product Note

Reacts with the heavy and light chains of rat IgG and the light chains of rat IgM. Pre-adsorbed with None. May react with immunoglobulins from other species and the light chains of other rat immunoglobulins

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pepsin digest of Rabbit Anti-Rat IgG(H+L).
Purification	Purified IgG
Conjugation	Tetramethylrhodamine (TRITC) Wavelength



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).