

Donkey Anti-Goat IgG antibody, pre-adsorbed (PE)

Cat. No. GTX04171-08

Host	Donkey
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Fr, FCM, ELISA, Multiplexing
Reactivity	Goat

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
ELISA	≤ 1 µg/mL
Multiplexing	Assay dependent

Not tested in other applications.

Product Note

Pre-adsorbed with Mouse and rat serum proteins (SP). May react with immunoglobulins from other species and the light chains of other goat immunoglobulins.

Properties

Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from donkeys hyperimmunized with goat IgG
Purification	Purified by antigen-affinity chromatography
Conjugation	Phycoerythrin (PE) Wavelength

For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).