

Mouse Anti-Chicken IgA antibody [A-1] (Biotin)

Cat. No. GTX04174-02

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ICC/IF, IHC-P, IHC-Fr, FCM, IP, ELISA
Reactivity	Chicken

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	$\leq 10 \mu\text{g/ml}$
IHC-P	Assay dependent
IHC-Fr	$\leq 10 \mu\text{g/ml}$
FCM	$\leq 1 \mu\text{g}/10^6\text{cells}$
IP	Assay dependent
ELISA	Assay dependent

Note : The suggested use of these reagents is in a final volume of 100µl.

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Affinity purified chicken Ig or isolated lymphocytes
Purification	Purified IgG2b
Conjugation	Biotin


For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX04174-02 IHC-Fr Image

IHC-Fr analysis of chicken intestine tissue using GTX04174-02 Mouse Anti-Chicken IgA antibody [A-1] (Biotin) followed by Streptavidin (CY3).



For full product information, images and publications, please visit our [website](#).