

Mouse Anti-Human IgG (Fc) antibody [H2] (PE)

Cat. No. GTX04177-08

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, FCM, ELISA, ELISPOT, Multiplexing
Reactivity	Human, Monkey, Chimpanzee

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	$\leq 0.1 \mu\text{g}/10^6\text{cells}$
ELISA	$\leq 1 \mu\text{g}/\text{mL}$
ELISPOT	Assay dependent
Multiplexing	Assay dependent

Note : The suggested use of these reagents is in a final volume of 100µl.

Not tested in other applications.

Product Note Reacts with Human/Rhesus/Cynomolgus/Chimpanzee/African Green Monkey/Mangabey/Pigtail Macaque IgG Fc.

Properties

Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified IgG2b
Conjugation	Phycoerythrin (PE) Wavelength




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx04177-08 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTx04177-08 Mouse Anti-Human IgG (Fc) antibody [H2] (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).