

## Mouse Anti-Human IgG1 (hinge region) antibody [4E3] (PE)

Cat. No. GTX04178-08

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, ELISA, Multiplexing, Purification
Reactivity	Human

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
FCM	Assay dependent
ELISA	≤ 1 µg/mL
Multiplexing	Assay dependent
Purification	Assay dependent

Not tested in other applications.

## Properties


Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified IgG1
Conjugation	Phycoerythrin (PE) <a href="#">Wavelength</a>

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTx04178-08 ELISA Image**

ELISA analysis of purified human immunoglobulins using serially diluted GTx04178-08 Mouse Anti-Human IgG1 (hinge region) antibody [4E3] (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).