

## Mouse Anti-Human IgG1 (Fc) antibody [HP6001] (AP)

Cat. No. GTX04179-03

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, FCM, ELISA, Multiplexing
Reactivity	Human

Package  
1 ml

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	Assay dependent
ELISA	1:2000-1:8000
Multiplexing	Assay dependent

Not tested in other applications.

## Properties

Form	Liquid
Buffer	50mM Tris, 1mM MgCl <sub>2</sub> , 50% Glycerol
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG1 myeloma protein
Purification	Purified IgG2b
Conjugation	Alkaline Phosphatase (AP)


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTx04179-03 ELISA Image**

ELISA analysis of purified human immunoglobulins using serially diluted GTx04179-03 Mouse Anti-Human IgG1 (Fc) antibody [HP6001] (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).