

## Mouse Anti-Human IgA2 (Fc) antibody [A9604D2] (PE)

**Cat. No. GTX04183-08**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	WB, IHC-P, IHC-Fr, FCM, ELISA, Depletion, ELISPOT, Multiplexing
<b>Reactivity</b>	Human

**Package**  
100 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
ELISA	≤ 1 µg/mL
Depletion	Assay dependent
ELISPOT	Assay dependent
Multiplexing	Assay dependent

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, a stabilizer
<b>Preservative</b>	0.09% sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
<b>Concentration</b>	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Purified human IgA2m(1) myeloma protein
<b>Purification</b>	Purified IgG1
<b>Conjugation</b>	Phycoerythrin (PE) <a href="#">Wavelength</a>




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

**DATA IMAGES**

**GTX04183-08 ELISA Image**

ELISA analysis of purified human immunoglobulins using serially diluted GTX04183-08 Mouse Anti-Human IgA2 (Fc) antibody [A9604D2] (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).