

Mouse Anti-Human IgE (Fc) antibody [B3102E8] (Biotin)

Cat. No. GTX04184-02

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-Fr, FCM, ELISA
Reactivity	Human

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
ELISA	1:5000-1:10000

Not tested in other applications.

Properties


Form	Liquid
Buffer	PBS
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified IgG1
Conjugation	Biotin

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04184-02 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTx04184-02 Mouse Anti-Human IgE (Fc) antibody [B3102E8] (Biotin) followed by Streptavidin (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).