

Mouse Anti-Human lambda light chain antibody [JDC-12] (HRP)

Cat. No. GTX04185-01

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, FCM, IP, ELISA
Reactivity	Human, Monkey

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	1:1000-1:2000

Not tested in other applications.

Product Note Reacts with Human/Rhesus Lambda.

Properties

Form	Liquid
Buffer	PBS, 50% glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Purification	Purified IgG1
Conjugation	Horseradish peroxidase(HRP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04185-01 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTX04185-01 Mouse Anti-Human lambda light chain antibody [JDC-12] (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).