

## Mouse Anti-Human lambda light chain antibody [JDC-12] (Biotin)

Cat. No. GTX04185-02

|              |                    |
|--------------|--------------------|
| Host         | Mouse              |
| Clonality    | Monoclonal         |
| Isotype      | IgG1               |
| Applications | WB, FCM, IP, ELISA |
| Reactivity   | Human, Monkey      |

Package  
500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution                  |
|--------------------|---------------------------------------|
| WB                 | Assay dependent                       |
| FCM                | $\leq 1 \mu\text{g}/10^6\text{cells}$ |
| IP                 | Assay dependent                       |
| ELISA              | 1:5000-1:20000                        |

**Note : The suggested use of these reagents is in a final volume of 100µl.**

Not tested in other applications.

**Product Note** Reacts with Human/Rhesus Lambda.

## Properties


|               |                                                                                         |
|---------------|-----------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                  |
| Buffer        | PBS                                                                                     |
| Preservative  | 0.09% sodium azide                                                                      |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| Concentration | 0.5 mg/ml (Please refer to the vial label for the specific concentration.)              |
| Purification  | Purified IgG1                                                                           |
| Conjugation   | Biotin                                                                                  |

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX04185-02 ELISA Image**

ELISA analysis of purified human immunoglobulins using serially diluted GTX04185-02 Mouse Anti-Human lambda light chain antibody [JDC-12] (Biotin) followed by Streptavidin (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).