

Mouse Anti-Human IgG4 (pFc') antibody [HP6023]

Cat. No. GTX04186

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Applications	WB, IHC-Fr, FCM, IP, ELISA
Reactivity	Human

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Properties


Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG4 myeloma protein
Purification	Purified IgG3
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04186 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTX04186 Mouse Anti-Human IgG4 (pFc') antibody [HP6023] followed by Goat Anti-Mouse IgG3 antibody, pre-absorbed (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).