

Mouse Anti-Human IgG4 (pFc') antibody [HP6023] (FITC)

Cat. No. GTX04186-06

| | |
|---------------------|----------------------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG3 |
| Applications | WB, IHC-Fr, FCM, IP, ELISA |
| Reactivity | Human |

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB | Assay dependent |
| IHC-Fr | Assay dependent |
| FCM | Assay dependent |
| IP | Assay dependent |
| ELISA | 1:200-1:400 |

Not tested in other applications.

Properties

| | |
|----------------------|---|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 0.09% sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light. |
| Concentration | 0.5 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | Human IgG4 myeloma protein |
| Purification | Purified IgG3 |
| Conjugation | Fluorescein isothiocyanate (FITC) Wavelength |


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04186-06 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTx04186-06 Mouse Anti-Human IgG4 (pFc') antibody [HP6023] (FITC).



For full product information, images and publications, please visit our [website](https://www.genetex.com).