

Mouse Anti-Human IgG4 (Fc) antibody [HP6025] (PE)

Cat. No. GTX04187-08

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, IHC-Fr, FCM, ELISA, Multiplexing, Purification
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
ELISA	≤ 1 µg/mL
Multiplexing	Assay dependent
Purification	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG4 myeloma protein
Purification	Purified IgG4 Fc
Conjugation	Phycoerythrin (PE) Wavelength


For full product information, images and publications, please visit our [website](#).

Date 2026 / 09 / 03 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX04187-08 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTX04187-08 Mouse Anti-Human IgG4 (Fc) antibody [HP6025] (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).