

Mouse Anti-Human IgG3 (hinge region) antibody [HP6050] (AP)

Cat. No. GTX04188-03

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, IHC-Fr, FCM, IP, ELISA, ELISPOT, Multiplexing, Purification
Reactivity	Human

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	1:500-1:1000
ELISPOT	Assay dependent
Multiplexing	Assay dependent
Purification	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	50mM Tris, 1mM MgCl ₂ , 50% Glycerol
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG3 myeloma protein
Purification	Purified IgG1
Conjugation	Alkaline Phosphatase (AP)




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx04188-03 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTx04188-03 Mouse Anti-Human IgG3 (hinge region) antibody [HP6050] (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).