

## Mouse Anti-Human IgG3 (hinge region) antibody [HP6050] (PE)

Cat. No. GTX04188-08

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, IHC-Fr, FCM, IP, ELISA, ELISPOT, Multiplexing, Purification
Reactivity	Human

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	≤ 1 µg/mL
ELISPOT	Assay dependent
Multiplexing	Assay dependent
Purification	Assay dependent

Not tested in other applications.

## Properties


Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG3 myeloma protein
Purification	Purified IgG1
Conjugation	Phycoerythrin (PE) <a href="#">Wavelength</a>

For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

**DATA IMAGES**

**GTX04188-08 ELISA Image**

ELISA analysis of purified human immunoglobulins using serially diluted GTX04188-08 Mouse Anti-Human IgG3 (hinge region) antibody [HP6050] (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).