

Goat Anti-Mouse IgG (Fab) antibody, pre-adsorbed

Cat. No. GTX04191

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Mouse

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	$\leq 1 \mu\text{g/mL}$

Not tested in other applications.

Product Note

Reacts with the Fab region of mouse IgG1, IgG2a, IgG2b, IgG2c, and IgG3. Pre-adsorbed with Mouse IgG Fc. May react with immunoglobulins from other species and the light chains of other mouse immunoglobulins

Properties

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from goat hyperimmunized with mouse IgG
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04191 ELISA Image

ELISA analysis of purified mouse immunoglobulins using serially diluted GTX04191 Goat Anti-Mouse IgG (Fab) antibody, pre-adsorbed followed by Mouse Anti-Goat IgG Fc antibody (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).