

## Goat Anti-Mouse IgM (Mu chain) antibody, pre-adsorbed (PE-Cy7)

**Cat. No. GTX04192-10**

<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, ICC/IF, IHC-P, IHC-Fr, FCM, IP, ELISA, Activation, EM, ELISPOT, Purification
<b>Reactivity</b>	Mouse

**Package**  
 250 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	$\leq 0.1 \mu\text{g}/10^6\text{cells}$
IP	Assay dependent
ELISA	Assay dependent
Activation	Assay dependent
EM	Assay dependent
ELISPOT	Assay dependent
Purification	Assay dependent

**Note : The suggested use of these reagents is in a final volume of 100µl.**

Not tested in other applications.

**Product Note**

Pre-adsorbed with Mouse IgG and IgA; human immunoglobulins and pooled sera. May react with immunoglobulins from other species.

## Properties


<b>Form</b>	Liquid
<b>Buffer</b>	PBS, a stabilizer
<b>Preservative</b>	0.09% sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.



For full product information, images and publications, please visit our [website](#).

<b>Concentration</b>	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Pooled antisera from goats hyperimmunized with mouse IgM
<b>Purification</b>	Purified by antigen-affinity chromatography
<b>Conjugation</b>	Phycoerythrin-Cyanine7 (PE-Cy7) <a href="#">Wavelength</a>
<b>Note</b>	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>

## DATA IMAGES



### GTx04192-10 FCM Image

FACS analysis of human peripheral blood lymphocytes using Mouse Anti-Human CD57 antibody followed by GTx04192-10 Goat Anti-Mouse IgM (Mu chain) antibody, pre-adsorbed (PE-CY7).



For full product information, images and publications, please visit our [website](#).